

YMC

NEW

LAUNCH



YMC HPLC & GC Vials
YMC Syringe Filters

**ECONOMICAL
RANGE**

**LOT TO LOT
CERTIFIED
RANGE**

**GC VIALS
RANGE**

**BOX TO BOX
CERTIFIED
RANGE**

**YMC
PRODUCT
RANGE**

ECONOMICAL RANGE



IN09SV-1:5CG:
2ml Flat Base Screw
Thread Clear Glass
Patch Vials




IN09SV-1:5AB:
2ml Flat Base Screw
Thread Amber Glass
Patch Vials



IN09PPS1: 9-425
Screw Cap, Red
PTFE / White Silicon
bonded pre-slit septa



IN09PPS2: 9-425
Screw Cap, Red PTFE/
White Silicon non
bonded non slit septa



IN09PPS3: 9-425
Screw Cap, Red PTFE/
White Silicon bonded
non slit septa



IN09PPS4: 9-425
Screw Cap, Red PTFE/
White Silicon non
bonded pre-slit septa

GC VIALS RANGE



IN20CG-2:0VFB:
20ml Flat Bottom,
Crimp Headspace,
Clear Glass Vial.



IN20CG-2:0RFB:
20ml Flat Bottom R2,
Crimp Headspace,
Clear Glass Vial.



IN20CG-2:0VRB:
20ml Round Bottom,
Crimp Headspace,
Clear Glass Vial.



IN20CG-2:0SRRB:
20ml Short Neck
Round Bottom R2,
Crimp Headspace,
Clear Glass Vial



IN20AC03SS:
20mm Tan PTFE/White
Silicon Septa with
Aluminium crimp cap



IN20AC03SW: 20mm
White PTFE/White
Silicon Septa with
Aluminium crimp cap



IN11AC03SS: 11mm
Red PTFE/White
Silicone Septa with
Aluminium crimp cap



IN11CG-2:0CV:
Clear Glass 11mm
Crimp Top Vial,
wide opening



IN11AB-2:0CV:
Amber Glass 11mm
Crimp Top Vial,
wide opening


LOT TO LOT CERTIFIED CLEAR GLASS




YPR-CG-CNC-02M: 2in 1 Kit: 1.5ml Short Thread Premium Vial, 12x32mm, clear glass, wide opening, with overwind-barrier ; 9mm UltraBond Combination Seal, bonded: Pre-slit



YPR-CG-CNC-02NS: 2in 1 Kit: 1.5ml Short Thread Premium Vial, 12x32mm, clear glass, wide opening, with overwind-barrier ; 9mm UltraBond Combination Seal, bonded : non slit.




YPR-CG-CNC-02PS: 2in 1 Kit: 1.5ml Short Thread Premium Vial, 12x32mm, clear glass, wide opening, with overwind-barrier ; 9mm Combination Seal, non bonded: pre slit.



YPR-CG-CNC-02NN: 2in 1 Kit 1.5ml Short Thread Premium Vial, 12x32mm, clear glass, wide opening, with overwind-barrier; 9mm Combination Seal, non bonded: non slit.




LOT TO LOT CERTIFIED AMBER GLASS




YPR-AB-CNC-02M: 2in 1 Kit:

1.5ml Short Thread Premium Vial,
12x32mm, amber glass, wide opening,
with overwind-barrier ; 9mm UltraBond
Combination Seal, bonded: slit.




YPR-AB-CNC-02NS: 2in 1 Kit:

1.5ml Short Thread Premium Vial,
12x32mm, amber glass, wide opening,
with overwind-barrier ; 9mm UltraBond
Combination Seal, bonded: non slit.



YPR-AB-CNC-02PS: 2in 1 Kit:

1.5ml Short Thread Premium Vial,
12x32mm, amber glass, wide opening,
with overwind-barrier ; 9mm
Combination Seal, non bonded: pre slit.



YPR-AB-CNC-02NN: 2in 1 Kit:


1.5ml Short Thread Premium Vial,
12x32mm, amber glass, wide opening,
with overwind-barrier ; 9mm
Combination Seal, non bonded: non slit.




BOX TO BOX CERTIFIED CLEAR GLASS




YPR-CG-CER-02M: LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread Premium Vial, 12x32mm, clear glass, wide opening, with overwind-barrier;9mm UltraBond Combination Seal, bonded:slit.



YPR-CG-CER-02NS: LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread Premium Vial, 12x32mm, clear glass, wide opening, with overwind-barrier;9mm Ultra Bond Combination Seal, bonded: non slit.




YPR-CG-CER-02PS: LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread Premium Vial, 12x32mm, clear glass, wide opening, with overwind-barrier;9mm Combination Seal, non bonded:pre slit




YPR-CG-CER-02NN: LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread Premium Vial, 12x32mm, clear glass, wide opening, with overwind-barrier;9mm Combination Seal, non bonded:non slit.




BOX TO BOX CERTIFIED AMBER GLASS




YPR-AB-CER-02M: LC/MS and GC/MS
Certified Vial Kit: 1.5ml Short Thread
Premium Vial, 12x32mm, amber glass,
wide opening, with overwind-barrier;9mm
UltraBond Combination Seal, bonded:slit.



YPR-AB-CER-02NS: LC/MS and GC/MS
Certified Vial Kit: 1.5ml Short Thread
Premium Vial, 12x32mm, amber glass,
wide opening, with overwind-barrier;9mm
Ultrabond Combination Seal, bonded:non slit.



YPR-AB-CER-02PS: LC/MS and GC/MS
Certified Vial Kit: 1.5ml Short Thread
Premium Vial, 12x32mm, amber glass,
wide opening, with overwind-barrier;9mm
Combination Seal, non bonded: pre slit.



YPR-AB-CER-02NN: LC/MS and GC/MS
Certified Vial Kit: 1.5ml Short Thread
Premium Vial, 12x32mm, amber glass,
wide opening, with overwind-barrier;9mm
Combination Seal, non bonded: non slit.





**NYLON
SYRINGE
FILTER**

Part No	Description	Media	Size	Pore Size	Pack Size / unit
INSFPA25045	YMC Nylon Syringe Filter, 13 mm, 0.2 μ m	Nylon	13 mm	0.2 μ m	100
INSFPA13045	YMC Nylon Syringe Filter, 13 mm, 0.45 μ m	Nylon	13 mm	0.45 μ m	100
INSFPA25020	YMC Nylon Syringe Filter, 25 mm, 0.2 μ m	Nylon	25 mm	0.2 μ m	100
INSFPA25045	YMC Nylon Syringe Filter, 25 mm, 0.45 μ m	Nylon	25 mm	0.45 μ m	100

**PTFE
SYRINGE
FILTER**

INSFPT13020	YMC PTFE Syringe Filter, 13 mm, 0.2 μ m	PTFE	13 mm	0.2 μ m	100
INSFPT13045	YMC PTFE Syringe Filter, 13 mm, 0.45 μ m	PTFE	13 mm	0.45 μ m	100
INSFPT25020	YMC PTFE Syringe Filter, 25 mm, 0.2 μ m	PTFE	25 mm	0.2 μ m	100
INSFPT25045	YMC PTFE Syringe Filter, 25 mm, 0.45 μ m	PTFE	25 mm	0.45 μ m	100

**PVDF
SYRINGE
FILTER**

INSFPD13020	YMC PVDF Syringe Filter, 13 mm, 0.2 μ m	PVDF	13 mm	0.2 μ m	100
INSFPD13045	YMC PVDF Syringe Filter, 13 mm, 0.45 μ m	PVDF	13 mm	0.45 μ m	100
INSFPD25020	YMC PVDF Syringe Filter, 25 mm, 0.2 μ m	PVDF	25 mm	0.2 μ m	100
INSFPD25045	YMC PVDF Syringe Filter, 25 mm, 0.45 μ m	PVDF	25 mm	0.45 μ m	100

Worldwide Availability

YMC America, Inc.
www.ymcamerica.com

YMC Europe GmbH
www.ymc.de

YMC Switzerland LLC
www.ymc-schweiz.ch

YMC Shanghai Rep. Office
www.ymcchina.com

YMC India Pvt. Ltd.
www.ymcindia.com

YMC Korea Co., Ltd.
www.ymckorea.com

YMC Taiwan Co., Ltd.
www.ymctaiwan.com

YMC Singapore Tradelinks Pte. Ltd.
www.ymc.co.jp/en/

YMC CO., LTD.

YMC Karasuma-Gojo Bld. 284 Daigo-cho,
Karasuma Nishiiru Gojo-dori Shimogyo-ku,
Kyoto 600-8106 Japan
Phone +81(0)75-342-4515, FAX +81(0)75-342-4550
www.ymc.co.jp

YMC (India) Private Limited

Head Office : A-154-155, 1st Floor, Eros Boulevard,
Hotel Crown Plaza, Mayur Vihar Phase-1,
New Delhi - 110091
Contact No : 011-45041601/701
E-mail : sales@ymcindia.com

YMC (India) Private Limited

Plot No : 78/A/6, Phase VI, Industrial Park,
Jeedimetala, Gajularamaram Village,
Quthbullapur, Medchal, Hyderabad.
Contact : +91 7893811309 / +91 7893811700
Email : labsupport@ymcindia.com